








# Surgical Guide for Breast Conserving Surgery



blue-dye injection for marking



## Products



## Features


- ① Guide the tumor excision range in surgery by internal marker (dye injection) and by external marker (drawing on the skin)
- ② Provide effective access route and excision guide through the segmentation of tumor margin from patient's MRI
- ③ Support a selective tumor range guide before and after the pre-operative neoadjuvant chemotherapy
- ④ Support precise tumor excision with no unnecessary margin cut, nor leaving unremoved malignant lesion

## Procedure

### 1. segmentation



### 2. 3D modeling



### 3. 3D printing final product

